

PRODUCT INFORMATION



EQUIPMENT MAINTENANCE

ERINOX TERMO

SPECIAL CLEANER



Range of use	Erinox Termo is a special cleaner for automated washing of endoscopes.		
Features	Erinox Termo product has been developed especially for automated washing of endoscopes. Erinox Termo can be used in most of automated washers of endoscopes.		
Product declaration	Anionic Surfactants < 5 %		washing effect
	Nonionic Surfactants < 5 %		washing effect
	2-Propanol < 5 %		isopropanol, solvent
	Sodiumcitrate < 5 %		dispersing agent
	Sodiumpolyacrylate < 5 %		dispersing agent
	Sodiumgluconate 5 - 15 %		dispersing agent
	Monopropyleneglycol 5 - 15 %		solvent
	pH of the dilution 9,5		
Instructions for use	Connect the Erinox Termo Drum with the automatic dosing system of the machine. Dosage 0.3 – 1 %. Preprocessing and washing according to the machine manufacturer's instructions. The endoscope will be connected to the washing and disinfecting cleaner. Choose the right washing and disinfecting program, machine dilutes both the detergent and disinfectant automatically. As a disinfectant use Erinox Termo or Thermoton Endo.		
Notice	Before using Erinox Termo make sure that the the product can be used in the type of disinfector washer of endoscopes you have.		
Package size / Product code	3 x 5 l canister	8141	EAN 64 1260 00 8141 1
	20 l canister	8151	EAN 64 1260 00 8151 0
Storage	Warm storage 3 years.		
Environmental impact	Packages are made of PE-HD plastic, which is environmentally friendly, recyclable and burnable material. The product is biological dissoluble.		
Country of origin	Finland		



08/08